VOCO**Rebilda LC****GB**

CE 0482

Instructions for use**Product description:**

Rebilda LC is a light-curing, radiopaque composite for core build-up with excellent mechanical properties.

Rebilda LC in the light blue shade shows the preparation margins.

Rebilda LC in the dentine shade can be used for translucent restorations.

Indications:

Core build-up on vital and non-vital teeth and blocking-out.

Method of use:**A) Drying**

Due to subsequent application of the light-curing core build-up material a dry tooth surface is required.

Any contamination of the cavity with blood or saliva needs to be prevented. The use of a rubber dam is recommended.

B) Take appropriate measures in proximity of pulp.**C) Bonding material**

Rebilda LC is used in the adhesive technique with a dentine/enamel bond.

Any light-curing bonding material may be used.

The product's instructions for use should be consulted to determine the correct preparation (etching technique) and procedure.

D) Application of Rebilda LC

Apply **Rebilda LC** with a suitable instrument.

It is recommended to use transparent matrices (smooth surfaces).

Place the tip of the polymerisation lamp as close as possible to the surface of the core build-up and polymerise with halogen light for at least 40 s.

If the distance is more than 5 mm the curing depth may be reduced.

Finishing of the core build-up can be conducted immediately after curing and removal of the shaping aids (e.g. with fine or extra-fine diamond burs, polishing discs) while using water-cooling.

Warnings/precautionary measures:

To achieve complete polymerisation avoid contact of **Rebilda LC** with phenolic substances (e. g. eugenol).

Insufficient curing of **Rebilda LC** may lead to pulp irritation.

Curing of larger core build-ups needs to be done layer by layer (approx. 2 - 4 mm each). When using light cones, transparent strips/matrices, extend polymerisation time correspondingly.

Insulate the core build-up thinly before taking impressions or constructing temporary crowns/bridges (it may be sufficient to moisten it well) to prevent adhesion between the materials.

Rebilda LC contains diurethane dimethacrylate, Bis-GMA, BHT, TEGDMA and barium aluminium borosilicate.

Rebilda LC should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Avoid contact with the oral mucosa.

The physical values of **Rebilda LC** correspond to DIN EN ISO 4049.

Information on storage and application:

Close **Rebilda LC** syringes after use.

Avoid exposure to heat. Store at temperatures of 4 °C to 23 °C.

Do not expose **Rebilda LC** to direct sunlight or the chairlight.

Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB1504 E1 0415 99 © by VOCO
VOCO GmbH