BISCO

CE 0459

CORE-FLO™ DC Lite

Dual-Cured

Low Viscosity Core Build-Up

Instructions for Use

BISCO, Inc.
1100 W. Irving Park Road
Schaumburg, IL 60193
U.S.A.
1-800-247-3368
1-847-534-6000
www.bisco.com

EC REP BISCO France
208, allée de la Coudoulette
13680 Lançon de Provence
France
33-4-90-42-92-92

IN-218R10
Rev. 5/24

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

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GENERAL INFORMATION

CORE-FLO DC Lite is a flowable dual-cured, radiopaque core build-up composite, specifically designed to cement posts and build-up the core structure in one easy application. **CORE-FLO DC Lite** is also an ideal dentin replacement material. **CORE-FLO DC Lite** in syringe and automix dual-syringe and #9 is available in both Natural/A1 & Opaque White shades. **CORE-FLO DC Lite** is intended for use with a bonding agent that is compatible with dual-cured materials, such as BISCO adhesives.

Intended user: licensed dentist.

Indications:
The principle indications for use of **CORE-FLO DC Lite** are:
1. Core build-up material
2. Cementation of posts or pins
3. Dentin replacement material

Warnings:
- Contamination by saliva, blood, or water should be avoided during these procedures. When using dental adhesives, contamination will substantially reduce bonding and may result in decreased longevity of the restoration.
- A rubber dam is recommended for all indications where contamination is possible.

Cautions:
Cross-contamination: Product may contain items that are designed for one-time use. Dispose of used or contaminated applicators, tips and mixers. Do not clean, sterilize or reuse.
- Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

Poisons:
- Avoid contact with the skin; unpolymerized (meth) acrylic resins may cause skin sensitization in some susceptible persons. In case of contact, wash skin with soap and water.
- To prevent polymerization or clogging, the dual-syringe mixing tip should not be stored with the unmixed pastes. If unfortunately, wipe mixing tip with disinfectant prior to storage.
- Unlike Bisco adhesives, do not use total-etch or self-etch systems bond to self-cured or dual-cured materials. Refer to the manufacturer's instructions.
- Isolate the cure build-up for taking impressions or fabricating temporary crown/veneer (it may be sufficient to moisten it thoroughly), in order to prevent any bonding with the materials.
- If refrigerated, product should be brought to room temperature before use.
- Protective eyewear should be worn by the patient and dental staff.
- Use only original component labels for specific expiration dates.
- Safety data sheet available at www.bisco.com
- Video: Use available upon request
- Instructions for Use available at www.bisco.com

TECHNICAL INFORMATION
Air room temperature (20°C/68°F - 25°C/77°F):
Working Time: Maximum 1 minute
Setting Time: Maximum 7 minutes
Working time and Setting time may vary based on storage conditions, temperature, humidity, etc.

INSTRUCTIONS FOR USE

CORE MATERIAL MIXING PROCEDURES:
1. Remove the syringe cap or mixing tip from the dual-syringe.
2. Dispense a small amount of material onto a mixing pad to check for even color and consistency of the materials.
3. Attach a mixing tip to the dual-syringe by lining up the tip with the syringe. Then turn the brown locking collar clockwise.
4. Pressing the plunger will mix and dispense the **CORE-FLO DC Lite** (See Section "DISPENSATION" Section for Working and Setting Time information.)

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