# 3M | Ketac™ Molar Easymix

ESPE

(en) **Glass Ionomer Restorative Material – Conventional**

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3M, ESPE, Ketac and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used in Germany. In Canada, © 2022, 3M. Unauthorized use prohibited. All rights reserved. Vita is a trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

MD CE 0123

3100019235/01

## (en) ENGLISH

**Product Description**
Ketac™ Molar Easymix is a radiopaque glass ionomer filling material for manual mixing and for use in bulk fillings. The material releases fluoride ions and can be used for luting purposes.
Ketac Molar Easymix is also suitable for use in minimal invasion filling therapy (MI) and for the sandwich technique. For the sandwich technique we recommend Ketac Molar Easymix is available in A3 shade based on the VITA Classical color system.
- These instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly identifiable. For safety and suitability of the mentioned products please refer to the corresponding Instructions for Use.

**Intended Purpose**
Conventional glass ionomer based dental restorative.
Intended users: educated dental professionals who have theoretical and practical knowledge on usage of dental products.
Clinical benefits: Restoration of initial and established tooth structure lesions.
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use.

**Indications**
Conventional filling therapy:
- As lining for single- or multi-surface composite fillings
- Core build-ups
- Primary tooth fillings
- Single-surface fillings in areas not involved in occlusion
- Class V fillings where appearance is not a priority
- Semi-permanent Class III and multi-surface fillings

Minimal intervention (MI) filling therapy (A.R.T. technique):
- Primary tooth fillings
- Single-surface fillings in areas not involved in occlusion
- Class V fillings where appearance is not a priority
- Semi-permanent Class III and multi-surface fillings
- Semi-permanent Class III fillings
- Fissure sealings

**Contraindications**
None

**Composition**
The ingredients of each component of the product are indicated in decreasing order of their concentration.
Powder: glass powder, polymeric acid. Liquid: water, polymeric acid, tartaric acid.
For more information, please consult the Safety Data Sheet/Safety Information Sheet (www.3M.com) or contact your 3M subsidiary.

**Precautionary Measures**
Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

**Potential Undesirable Side Effects and Complications**
Potential effects are local irritations. Restoration debonding or fracture can lead to swallowing or aspiration of the restoration or parts of it in rare cases.
The summary of safety and clinical performance for the device can be retrieved in the European database on medical devices (Eudamed) under the following link (as soon as the corresponding database module is available): https://ec.europa.eu/tools/eudamed
The Instructions for Use and the summary of safety and clinical performance are also available at: https://hcbgregulatory.3m.com
Basic-UDI-DI: 08082232761020000004EE2

**Steps Prior to Preparation**
Conventional filling therapy and A.R.T. technique:
- Keep handy all required instruments and materials.
- In the case of multi-surface fillings, place a matrix band to protect neighboring teeth and simplify contour shaping.

Especially in the A.R.T. technique:
A minimal set of instruments should include: dental mirror, probe, tweezers, excavators (large, medium, small), modeling instrument, carver (axe-shaped cutting instrument), plastic spatula, mixing pad, a matrix band or crown band or crown rolls, cotton pellets, water cup (sterilizable or disposable). This basic set can be supplemented with additional instruments depending on the available financial resources and local circumstances.
- Prior to commencing preparation, the working area should be thoroughly cleaned and dried. This renders discolored, demineralized enamel easy to recognize and treat.
- If no suction is available during treatment, the working area must be kept dry with repeated changes of cotton rolls.

**Cavity Preparation**
Conventional filling therapy:
- Remove only the carious substance; no undercuts required.
  - Do not prepare thin edges since the material requires a minimal wall thickness of 0.5 mm!
- Clean the cavity with water and dry with air.

A.R.T. technique:
- Extend the orifice of small cavities with a dental hatchet or carver to gain better access. As a consequence, any enamel lying soft will break off and a sufficient space for access of the smallest manual excavator is created.
- Commencing at the enamel-dentin interface, use circular abrasive motions of the smallest excavator to completely remove any carious dentin from the cavity. Unsupported enamel should also be removed.
- Subsequently clear the floor of the cavity.
  - Use the largest suitable excavator to treat large cavities, and exert slight pressure in order not to inadvertently expose the pulp due to excessive pressure!
- Clean the cavity with water and dry.

**Protection of the Pulp**
Conventional filling therapy and A.R.T. technique:
If necessary, cover the areas in close proximity to the pulp by applying spots of a hard-setting calcium hydroxide material to avoid any pulpal irritation.

**Conditioning**
Conventional filling therapy and A.R.T. technique:
- For optimal chemical bonding to the teeth, the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner, or, for use with the A.R.T. technique, Ketac Molar Easymix Liquid: to the prepared surfaces and allowing the substance(s) to react for 10 sec.
- Then rinse with copious amounts of water and blow dry in only 2–3 short intervals with water-free and oil-free air or dab dry with cotton pellets.
  - **Your cavity should have a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid further contamination; repeat the procedure if surfaces become recontaminated.

**Times**
At an ambient temperature of 23 °C/73 °F and 50% relative humidity, the following times apply:

| | min:sec |
|---|---|
| Mixing | 00:30 |
| Processing from start of mixing | 02:00 |
| Setting from start of mixing | 05:00 |

At higher temperatures, the time available for processing is shortened while it is prolonged at lower temperatures (e.g., mixing on a chilled glass slab). The processing time is also shortened if the amount of powder exceeds the recommendation. Exceeding the indicated processing time may lead to reduced bonding to enamel and dentin.

**Dosing**
The mixing ratio by weight is 4.5 parts powder (1 level spoonful) to 1 part liquid (1 drop).
Conventional filling therapy and A.R.T. technique:
- Loosen the powder by shaking the bottle.
- Dip the enclosed spoon into the powder and use the plastic scraper at the bottle opening to level it off. Avoid compressing the powder.
- Dose a sufficient amount of powder on a mixing pad.
- Hold the bottle of liquid in a vertical position and dose the appropriate number of drops for the amount of powder next to the powder.
  - There must not be any dried residue of the liquid on the nozzle, and the bottle must be closed properly after use!
- Carefully re-close both bottles after use and ensure that no powder residue remains on the plastic scraper of the powder bottle.

**Mixing**
Conventional filling therapy and A.R.T. technique:
- Process Ketac Molar Easymix at room temperature (20–25 °C/68–77 °F).
- Use a metal or plastic cement spatula and a mixing pad or glass slab for mixing.
- Ideally, the powder should be transferred into the liquid in not more than two portions.
- Break out the paste repeatedly until the consistency is homogeneous.
  - Moisture can have a negative effect on the properties of the glass ionomer cement. Therefore, ensure pastes is not obtained with the given mixing ratios, discard the material.

**Application**
Avoid any and all contamination from water, saliva, etc. and keep the working area dry during the entire time of the application and shaping.
Conventional filling therapy:
- Apply Ketac Molar Easymix in several portions using a suitable instrument. Undercuts may need to be treated first.
  - Carefully avoid trapping air bubbles under enamel undercuts.

A.R.T. technique:
- Apply Ketac Molar Easymix in several portions using a suitable instrument. Undercuts may need to be treated first.
  - Carefully avoid trapping air bubbles under enamel undercuts.
  - Slightly overfill cavities; apply the material to adjacent grooves and fissures.
- Apply a small amount of Vaseline to your forefinger and use this finger to firmly press the glass ionomer material into cavities, fissures, and grooves.
- Repeatedly move the tip of the forefinger in a mesio-distal and bucco-lingual direction to distribute the material over the entire occlusal surface (finger pistol).
- Finally, remove the finger sideways to prevent the restorative material from lifting out of the cavity.

Fissure sealing with the A.R.T. technique:
The same steps of cleaning, conditioning, application, and processing as in the A.R.T. technique are performed to seal fissures, the only difference being that there is no preparation involved.

**Filling Protection**
Conventional filling therapy and A.R.T. technique:
The application of a protective layer on the filling is not necessary. In the event that a protective layer is desired, varnishes and/or cotton rolls can be used after filling protection as described below:
- Apply vaseline to all of the exposed surfaces of the glass ionomer cement immediately after shaping.

**Finishing**
Conventional filling therapy:
- Beginning no earlier than 5 minutes after the start of mixing, use Arkansas stone, fine-grain diamonds, polishing discs of decreasing grain size, or silicone polishing cups and points.
- If desired, apply vaseline again for surface sealing.
- The patient must not expose the filling to any pressure for 1 hour!

A.R.T. technique:
- Remove any excess material with a carver or large excavator as soon as possible without pulling the material from the cavity.
- Before the filling material is completely cured, test the occlusion with some articulation paper.
- Carefully remove any increase in vertical dimension with a carver.
- Especially at the edges of the restoration, the surface should be slightly out-of-contact.
- Apply more Vaseline to the working areas to protect the restoration.

The patient must not expose the filling to any pressure for 1 hour!

**Cleaning and Disinfection**
**Caution**
Make sure that the cleaning and disinfecting agents you have chosen do not contain any of the following materials:
- Oxidizing agents (e.g., hydrogen peroxides)
- Alkalis
- Glutaraldehydes

Use only validated methods for cleaning and disinfection.
Use solely cleaning and disinfecting agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). CavWipes™ have been tested for use as cleaning and disinfecting agent. Always observe all applicable legal and hygiene regulations for control offices and/or hospitals.

**General**
The following procedure should be applied after use or before any subsequent use. The glass bottle/plastic bottle/dosing spoon must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

**Manual Cleaning by Wiping**
- Clean the glass bottle/plastic bottle/dosing spoon after or before every application with ready-to-use disinfectant wipes (e.g., CavWipes™ based on the active substances alcohols and quaternary compounds) or equivalent cleaning agents. Follow the cleaning agent manufacturer's directions for use.
- Thoroughly clean the surface of the glass bottle/plastic bottle with a new ready-to-use CavWipe™ and the dosing spoon with two ready-to-use CavWipes™ for at least 30 seconds or until no more dirt is visible.
- Special care must be exercised in cleaning of the grooves on the lid and recess at the tip.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**
- Disinfect the bottle with ready-to-use disinfectant wipes (e.g., CavWipes™ based on the active substances alcohols and quaternary compounds) or equivalent disinfection agents. Follow the disinfecting agent manufacturer's directions for use, observe the exposure time.
- Disinfect the surface of the glass bottle for three minutes using at least two new disinfectant wipes. Keep all surfaces moist during the six-minute period. Special care must be exercised in disinfection of the grooves on the lid.
- Disinfect the surface of the dosing spoon three minutes using at least one new disinfectant wipe. Keep all surfaces moist during the three-minute period. Special care must be exercised in disinfection of grooves or recess at the logo.
- Note: Cleaning residues on the dosing spoon: rinse the spoon at least 10 sec. with tap water and dry with a clean, soft, lint free cloth.

**Inspection, Maintenance, Testing**
- Examine the bottles/dosing spoon for damage, discoloration, and contamination after each cleaning.
- Do not use damaged devices under any circumstances. If the bottles/dosing spoon are not visibly clean, repeat the cleaning and disinfection procedure.

**Notes**
- The material bonds to metal instruments. Residual quantities of the cement should be removed with water before they set on the instrument.
- Ketac Molar Easymix is a multiple patient multiple use device.

**Disposal**
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of the substance(s) to react for 10 sec.

**Storage and Shelflife**
Store the product at 2–25 °C/59–77 °F.
**To protect the powder from moisture, do not open the powder bag until immediately prior to the first application.**
**Replace the cap securely after use.**
Do not use after the expiration date.

**Customer Information**
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

**Limitation of Liability**
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

**Symbol Glossary**

| Reference Number and Symbol Title | Symbol | Description of Symbol |
|---|---|---|
| ISO 15223-1 5.1.1 Manufacturer | | Indicates the medical device manufacturer. |
| Swiss Authorized Representative | CH REP | Indicates the authorised representative in Switzerland. |
| ISO 15223-1 5.1.3 Date of Manufacture | | Indicates the date when the medical device was manufactured. |
| ISO 15223-1 5.1.4 Use-by date | | Indicates the date after which the medical device is not to be used. |
| ISO 15223-1 5.1.5 Batch code | LOT | Indicates the manufacturer's batch code so that the batch or lot can be identified. |
| ISO 15223-1 5.1.6 Catalogue number | REF | Indicates the manufacturer's catalogue number so that the medical device can be identified. |
| ISO 15223-1 5.3.7 Temperature limit | | Indicates the temperature limits to which the product can be safely exposed. |
| ISO 15223-1 5.4.3 Consult instructions for use | | Indicates the need for the user to consult the instructions for use. |
| CE Mark | CE | Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement. |
| CE Mark | CE | Indicates conformity to all applicable European Union Medical Device Regulations and Directives. |
| ISO 15223-1 5.7.7 Medical Device | MD | Indicates the item is a medical device. |
| Corrosive | | Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dentist, physician, professional, Or Code of Federal Regulations (CFR) sec. 801.109(b)(1) |
| Low-density polyethylene | LDPE | Indicates product packaging is made of low-density polyethylene. Official Journal of the EC, Commission Decision (97/129/EC) |
| Aluminum and Low-density polyethylene | ALU LDPE | Indicates product packaging is made of aluminum and low-density polyethylene. Official Journal of the EC, Commission Decision (97/129/EC) |
| Colorless Glass | GL | Indicates product packaging is made of colorless glass. Official Journal of the EC, Commission Decision (97/129/EC) |
| Green Dot | | Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe. |
| ISO 15223-1 5.1.8 Importer | | Indicates the entity importing the medical device into the locale. |

For more information, see HCBGregulatory.3M.com
Information valid as of August 2022

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3M, ESPE, Ketac and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.