VOCO

Registrado® X-tra Clear

EN Instructions for use

MD EU Medical device

CE 0482

Product description:
Registrado X-tra Clear is an addition-curing silicone impression material, particularly for the composite injection molding technique, which stands out for its transparency, high final hardness of 70 Shore A, and a fast setting time.

Indications:
- Silicone matrix for the composite injection molding technique
- Bite registration
- Occlusal surface stamps for light-cured composite registrations
- Tool for fabricating temporary veneers
- Transfer matrix for attaching brackets

Contraindications:
Registrado X-tra Clear contains platinum catalysts, peppermint flavor and methylpentynol. **Registrado X-tra Clear** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
Registrado X-tra Clear is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
Registrado X-tra Clear should only be used by a professionally trained dental practitioner.

Use:
Usable mixing tips: type 16

1. Preparation:
To dispense **Registrado X-tra Clear**, use a **VOCO** Dispenser type 2 or any other suitable dispenser. For technical reasons there may be slight differences in the fill levels of both chambers. You should therefore remove the cap and press out the material until an equal amount has been ejected from both outlets. Then attach the mixing tip and lock (with a 90° turn clockwise). The material is automatically mixed by the mixing tip in the correct ratio.
The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories. The attached mixing tip is intended for single use only. After use, store the cartridge tightly sealed. Check that material flows freely through the openings before reusing.

2. Processing:
To obtain an occlusal registration, **Registrado X-tra Clear** is applied directly to the patient's upper or lower dental arch. Bite down within 30 seconds. At room temperature (23°C), **Registrado X-tra Clear** has a processing time of max. 45 seconds. Following an intraoral retention time of 2 minutes after closing the mouth, the occlusal registration can be removed. When using extraorally, please note the longer curing time of 5 min. After curing, **Registrado X-tra Clear** can be finished in a simple and controlled manner using a scalpel or scissors. It remains dimensionally stable in the long term.

Technical data:

Processing time (at 23°C)	≤ 45 sec
Retention time in mouth (at 37°C)	≥ 2 min
Minimum dwell time on a model (at 23 °C)	≥ 5 min
Shore A hardness	approx. 70

Warnings, precautionary measures:
- Do not use **Registrado X-tra Clear** in combination with condensation-curing silicones.
- The setting behaviour of **Registrado X-tra Clear** can be affected by the use of latex gloves, hand creams, cleaning agents, mouth rinses, retraction solutions and other material residue (e.g., impression materials). Thus, these must be removed prior to application (spray clean if necessary).
- Avoid contact between incompletely cured silicone material and clothing.
- Disinfection is recommended prior to further processing of the occlusal registration created.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Disinfection:
Bite registrations can be treated with a 2 % glutaraldehyde solution, a 5.25 % sodium hypochlorite solution, or another disinfectant solution suitable for A-silicones. The instructions for use provided by the respective disinfectant manufacturer must be followed.

Storage:
Storage at 4°C-23°C. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Registrado X-tra Clear** must be reported to VOCO GmbH and the responsible authority.

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VC 60 AA2030 E1 0726 99 © by VOCO